



BrdU incorporation in multiparameter flow cytometry : A new cell cycle assessment approach in multiple myeloma

Journal:	<i>Cytometry: Part B - Clinical Cytometry</i>
Manuscript ID	18-031.R2
Wiley - Manuscript type:	Original Article
Date Submitted by the Author:	n/a
Complete List of Authors:	Requirand, Guilhem; Centre Hospitalier Regional Universitaire de Montpellier Robert, Nicolas; Centre Hospitalier Regional Universitaire de Montpellier Boireau, Stephanie; Centre Hospitalier Regional Universitaire de Montpellier Vincent, Laure; Centre Hospitalier Regional Universitaire de Montpellier Seckinger, Anja; University of Heidelberg, Innere Medizin V Bouhya, Salaheddine; Centre Hospitalier Regional Universitaire de Montpellier Ceballos, Patrice; Centre Hospitalier Regional Universitaire de Montpellier Cartron, Guillaume; Centre Hospitalier Regional Universitaire de Montpellier Hose, Dirk; University of Heidelberg, Innere Medizin V Klein, Bernard; Université Montpellier, CHU Montpellier, Institute of Research in Biotherapy Moreaux, Jerome; University of Montpellier, UFR Medecine
Key Words:	Multiple Myeloma
Controlled Key Words:	proliferation, prognosis, Multiparameter analysis, Cytometry, Monoclonal

SCHOLARONE™
Manuscripts

1
2
3 **BrdU incorporation in multiparameter flow cytometry : A new cell cycle**
4 **assessment approach in multiple myeloma**

5
6 Guilhem Requirand^{1,2}, Nicolas Robert^{1,2}, Stéphanie Boireau^{1,2}, Laure Vincent⁶, Anja
7 Seckinger^{4,5}, Salaheddine Bouhya⁶, Patrice Ceballos⁶, Guillaume Cartron^{3,6}, Dirk
8 Hose^{4,5}, Bernard Klein^{1,2,3}, Jérôme Moreaux^{1,2,3}
9

10
11 ¹ CHU Montpellier, Department of Biological Hematology, Montpellier, France

12 ² IGH, CNRS, Univ Montpellier, France

13 ³ Univ Montpellier, UFR de Médecine, Montpellier, France

14 ⁴ Medizinische Klinik und Poliklinik V, Universitätsklinikum Heidelberg, Heidelberg,
15 Germany

16 ⁵ Nationales Centrum für Tumorerkrankungen, Heidelberg, Germany

17 ⁶ CHU Montpellier, Department of Clinical Hematology, Montpellier, France.

18 ⁷ Univ Montpellier, UMR CNRS 5235, Montpellier, France.
19
20
21
22
23
24

25 Corresponding author: Jerome MOREAUX,

26 Laboratory for Monitoring Innovative Therapies

27 Department of Biological Hematology

28 Hôpital Saint-Eloi - CHRU de Montpellier

29 80, av. Augustin Fliche

30 34295 Montpellier Cedex 5

31 IGH - Institute of Human Genetics

32 CNRS-UM UMR 9002

33 <http://www.igh.cnrs.fr/>

34 phone: +33 (0)4 67 33 79 03

35 fax: +33 (0)4 67 33 70 36

36 mail: jerome.moreaux@igh.cnrs.fr
37
38
39
40
41
42
43

44 **Conflict of interest:**

45 The authors have no conflict of interest to declare.
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

Abstract

Background

Multiple myeloma (MM) is a neoplasia characterized by the accumulation of malignant plasma cells (PC) in the bone marrow. Although proliferation markers have been studied in MM, none of the current staging systems include them. Moreover, approaches used to analyze proliferation do not separate MM cells (MMCs) from normal PC.

Methods

In this study, we combined multiparameter flow cytometry and BrdU incorporation or Ki67 staining to analyze MM cell proliferation in 44 monoclonal gammopathy of undetermined significance (MGUS), 153 newly diagnosed MM patients and 69 MM patients at relapse. The prognostic value of proliferation assessment was analyzed in 60 newly diagnosed patients treated with high-dose chemotherapy supported by autologous hematopoietic stem cell transplantation.

Results

The median number of proliferating malignant PC significantly increases during MM disease progression. MM patients with a percentage of proliferating MMCs > 1.42% using BrdU/DAPI or > 1.1% using ki67/DAPI, are associated with a significantly shorter event free survival compared to patients with a lower percentage of proliferating MMCs.

Conclusions

Combination of flow cytometry with BrdU or ki67/DAPI staining could become a standard for the determination of MM cell proliferation. Furthermore, in the context of new effective myeloma treatment options, assessment of MM cell proliferation may be valuable, in clinical trials, to identify novel agents that could significantly affect the

1
2
3 small proliferative compartment of MM cells.
4
5
6
7
8
9

10 **Introduction**

11
12
13 Multiple myeloma (MM) is characterized by an accumulation of malignant plasma
14 cells (PC) in the bone marrow (BM) leading to bone destruction, renal failure and
15 organ dysfunction. Within the last decade, new treatment strategies have significantly
16 improved the outcome of MM patients (1). However, despite these new treatments
17 using drugs such as proteasome inhibitors or immunomodulatory agents (IMiDs),
18 most patients will eventually relapse and die (2). A key biological question is the
19 ability to perform an accurate appreciation of the tumor burden, in terms of infiltrating
20 plasma cells, for the initial staging of patients and for the monitoring of the disease.
21 Flow cytometry is widely used for the diagnosis and monitoring of hematological
22 disorders including MM (3-5). Proliferating PC, i.e. the growing fraction of MM cells,
23 have been evaluated by flow cytometry detection of newly synthesized DNA using
24 techniques based on the uptake of tritiated thymidine or bromodeoxyuridine (BrdU)
25 (6-8). The so-called plasma cell labeling index (PCLI) based on propidium iodide
26 incorporation (9) has been shown to be a powerful and independent predictor of
27 survival in MM (10,11) as the assessment of Ki-67 expressing cells in the malignant
28 fraction (12). Although proliferation markers have been studied in MM, none of the
29 current prognostic scores or staging systems, available to physicians nowadays,
30 include these markers. Moreover, the approaches used to analyze proliferation in
31 MM do not separate MMCs from their normal cell counterparts. Therefore, easy-to-
32 perform markers of proliferation remain to be developed for routine practice. Novel
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

1
2
3 multiparameter flow cytometry (MFC) enables increased sensibility and depth of
4 malignant PC detection (13). Furthermore, assessment of proliferating myeloma cells
5 is of special interest, as proliferating MMCs can be targeted by available treatments
6 and upcoming therapeutic treatment options (14). The study aimed to combine MFC
7 and BrdU incorporation or Ki67 staining to analyze MM cell proliferation in samples
8 from MGUS and MM patients. We also investigated the prognostic value of
9 proliferation assessment in a cohort of newly diagnosed MM patients.
10
11
12
13
14
15
16
17

18 **Materials and methods**

19 **Primary multiple myeloma cells**

20
21
22
23 Bone marrow samples were collected after patients' written informed consent in
24 accordance with the Declaration of Helsinki and institutional research board approval
25 from Montpellier University hospital (DC-2008-417). In particular, bone marrow
26 samples were collected from 44 MGUS patients, 153 MM patients at diagnosis and
27 69 at relapse. For 60 newly diagnosed patients treated with high dose Melphalan
28 (HDM) and autologous stem cell transplantation (ASCT), MMCs were purified using
29 anti-CD138 MACS microbeads (Miltenyi Biotech, Bergisch Gladbach, Germany) and
30 their gene expression profile (GEP) obtained using Affymetrix U133 plus 2.0
31 microarrays.
32
33
34
35
36
37
38
39
40
41
42
43

44 **Antibodies**

45
46
47 Anti CD45 R Phycoerythrin-TexasRed (ECD) (clone-J33), anti CD138
48 allophycocyanin (APC) (clone-BB4), anti CD38 Pe-Cyanin 5.5 (PeCy5.5) (clone-
49 LS198-4-3), CD20 allophycocyanin-AlexaFluor 750 (APC-AF 750) (clone-B9E9),
50 CD19 R Phycoerythrin-Cyanine 7 (PC7) (J3-119), CD56 R Phycoerythrin-Cyanine 7
51 (PC7) (clone-N901), CD117 R Phycoerythrin-Cyanine 7 (PC7) (clone-104D2D1),
52
53
54
55
56
57
58
59
60

1
2
3 CD27 R Phycoerythrin-Cyanine 7 (PC7) (clone 1A4CD27) monoclonal antibodies
4 (mAbs), were purchased from Beckman Coulter (Miami, FL, USA), CD200
5
6 Phycoerythrin-Cyanine 7 (clone-MRCOX-104), Kappa immunoglobulin light chain
7
8 phycoerythrin (PE) (clone-G20-193), Lambda immunoglobulin light chain fluorescein
9
10 isothiocyanate (FITC) (clone-JDC-12), Ki67 allophycocyanin (APC) (clone-B56) and
11
12 anti BrdU-APC (BrdU Flow Kit from BD Biosciences) mAbs were purchased from
13
14 Becton Dickinson Biosciences (San Jose, CA, USA).
15
16

17 18 **Cell staining for immunophenotyping of MM**

19
20
21 BM aspirates were harvested in heparin-containing syringes. Samples were filtered
22
23 through a 100 µm filter to remove bone debris. Whole bone marrow was incubated
24
25 with or without BrdU (10µM), (Flow Kit APC, Biosciences, San Jose, CA, USA) at
26
27 37°C during two hours. To lyse red blood cells, samples were incubated with a 4-fold
28
29 volume excess of ammonium chloride for 20 minutes. Then, leukocytes were washed
30
31 and incubated in phosphate buffered saline (PBS) (Sigma-Aldrich, Dorset, United
32
33 Kingdom) containing 10% goat serum albumin (PAA Laboratories, Paris, France) for
34
35 10 minutes as blocking step. Erythrocytes-lysed bone marrow samples after
36
37 saturation step using goat serum (PPA laboratories, Austria), were labeled with
38
39 optimal concentrations of antibodies (Table 1). Cells were stained with anti CD19,
40
41 anti CD20, anti CD27, anti CD38, anti CD56, anti CD200 and anti CD45 mAbs,
42
43 during twenty minutes at 4°C, washed twice in PBS/10% goat serum, fixed and
44
45 permeabilized (BrdU Flow Kit from BD Biosciences, BD Biosciences) and washed
46
47 using PermWash buffer (BD Biosciences). Cells were incubated during one hour at
48
49 37°C in PBS containing DNase I (300 µg/ml). Cells were then stained with anti-
50
51 Kappa light chain, anti-Lambda light chain and anti BrdU-APC mAb (BrdU Flow Kit
52
53 from BD Biosciences) during twenty minutes at 4°C (Table 1). For Ki67 staining, cells
54
55
56
57
58
59
60

1
2
3 were labeled as described in Table 1. After fixation and permeabilization, cells were
4
5 stained with an APC conjugated anti-KI67 antibody or isotypic APC (BD Biosciences)
6
7 and anti-Kappa light chain, anti-Lambda light chain mAb. Nuclei were stained in
8
9 PermWash containing DAPI (2 μ g/ml) (4',6 Diamido-2-phenylindole, dilactate from
10
11 Invitrogen, Life Technologies) during twenty minutes at 4°C. Cells were analyzed
12
13 using Cyan cytometer (Beckman Coulter, FL, USA). Contaminating events were
14
15 removed on both FCS and SSC plots. Singlets were plotted on FSC-A vs FSC-H and
16
17 SSC-A vs SSC-H plots to remove debris and to select the total leukocyte population.
18
19 PC and B cells were selected on CD45/CD38 and CD20/CD19 plots, respectively.
20
21 Abnormal PC were selected based on the CD19, CD27, CD56, CD117 and CD200
22
23 signals (13). 5,000,000 events were acquired per tube and the minimum number of
24
25 abnormal plasma cells needed was 20. The maximum sensitivity of the method is
26
27 0.0004% (13). The analyses have been done using FlowJo software (FlowJo, LLC).
28
29
30

31 **Gene expression profiling and statistical analyses**

32
33
34 Gene expression data were normalized with the MAS5 algorithm and analyses
35
36 processed with GenomicScape (<http://www.genomicscape.com>) (15) the R.2.10.1
37
38 and bioconductor version 2.5 programs (16). Gene Set Expression Analysis (GSEA)
39
40 was used to identify genes and pathways differentially expressed between
41
42 populations (17). Univariate and multivariate analysis of genes prognostic for
43
44 patients' survival was performed using the Cox proportional hazard model. Difference
45
46 in event free survival between groups of patients was assayed using Maxstat
47
48 algorithm (18) and survival curves plotted using the Kaplan-Meier method.
49
50
51

52 **Results**

1
2
3 MFC was used to delineate malignant and normal plasma cells (13) (Supplementary
4 Figure S1A). PC were selected on CD45/CD38 plot. Specific malignant markers
5 including CD117, CD56 and CD200 combined with CD19, CD27 normal PC markers
6 and kappa or lambda light chains were used to identify malignant from normal
7 plasma cells (Supplementary Figure 1A) (Table 1). In order to assess proliferation of
8 myeloma cells specifically, we combined MFC with BrdU incorporation or Ki67
9 staining (Table 1) (Supplementary Figure S1B).

10
11 MFC and BrdU incorporation was investigated in 44 MGUS, 153 newly diagnosed
12 MM patients and 69 MM patients at relapse (Figure 1A). The median percentage of
13 malignant PC in the S phase of the cell cycle in MGUS, newly diagnosed MM
14 patients and relapsing patients was 0.2% (range: 0 to 1.23%), 0.4% (range: 0 to
15 17.3%) and 0.83% (range: 0 to 33.8%) respectively. The median number of
16 proliferating malignant PC was 2-fold higher in newly diagnosed patients compared
17 to MGUS ($P = 3E-5$) and 2-fold higher in MM at relapse compared to patients at
18 diagnosis ($P = 0.03$) (Figure 1A). We also investigated the interest of Ki67 staining to
19 assess proliferation of malignant plasma cells. MFC incorporating Ki67/DAPI staining
20 was also investigated in 11 MGUS, 49 newly diagnosed MM patients and 18 MM
21 patients at relapse (Figure 2A). The median percentage of Ki67⁺/DAPI⁺ (cells in S-
22 G2-M phases) malignant PC in MGUS, newly diagnosed MM patients and relapsing
23 patients was 0.38% (range: 0.07 to 1.75%), 0.59% (range: 0 to 18.8%) and 0.91%
24 (range: 0 to 17.6%) respectively. The median number of Ki67⁺/DAPI⁺ malignant PC
25 was significantly higher in newly diagnosed patients compared to MGUS and in MM
26 at relapse compared to MGUS ($P = 0.04$). MMC proliferation assessment using BrdU
27 incorporation or Ki67⁺/DAPI⁺ staining was highly correlated ($r = 0.95$, $P < 0.001$)
28 (Supplementary Figure S2).

1
2
3 The prognostic value of proliferation assessment using BrdU incorporation was
4 analyzed in 60 newly diagnosed patients treated with high dose chemotherapy
5 supported by autologous hematopoietic stem cell transplantation (HSCT). Using R
6 Maxstat algorithm(18), allowing to determine the optimal cutpoint for continuous
7 variables (19-21), we identified that MM patients with a percentage of proliferating
8 MMCs > 1.42 are associated with a significant shorter event free survival (20 months
9 median EFS) ($P = 0.0005$) compared to patients with a percentage of proliferating
10 MMCs ≤ 1.42 (not reached median EFS) (Figure 1B). Using the same methodology,
11 MM patients with a percentage of ki67+/DAPI+ MM cells > 1.1% were significantly
12 associated with an adverse EFS (19 months median EFS) compared to patients with
13 a percentage of ki67+/DAPI+ MM cells $\leq 1.1\%$ (not reached median EFS) ($P = 3E-5$)
14 (Figure 2B). Therefore, proliferation staining using BrdU incorporation or ki67/DAPI
15 was significantly higher in newly diagnosed patients classified in the proliferation
16 subgroup of the myeloma molecular classification(22) (Supplementary Figure S3).
17 MM cell proliferation using BrdU incorporation or KI67/DAPI was also significantly
18 higher in patients with high growth proliferation gene index (GPI)(23), UAMS-high risk
19 score (HRS)(24), high RS GEP-based score(25) and high IFM score(26) associated
20 with an adverse outcome in MM (Supplementary Figure S4 and Supplementary
21 Figure S5).

22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45 Gene Set Expression Analysis (GSEA) was performed comparing gene expression
46 profiles of patients with high risk proliferation assessment (S phase > 1.42%) to
47 patients with low risk proliferation (S phase < 1.42%). Accordingly, high risk patients
48 are characterized by a significant enrichment of genes involved in cell cycle and
49 proliferation (CHANG CYCLING GENES, MOLENAAR TARGETS OF CCND1 CDK4,
50 ROSTY CERVICAL CANCER PROLIFERATION CLUSTER and CHIANG LIVER
51
52
53
54
55
56
57
58
59
60

1
2
3 CANCER SUBCLASS PROLIFERATION UP, $P < 0.005$) (Figure 1C). In contrast,
4
5 MM patients with low risk proliferation assessment were characterized by a
6
7 significant enrichment of genes related to mature bone marrow plasma cells
8
9 (MOREAUX B LYMPHOCYTE MATURATION BY TACI UP, $P = 0.002$) (Figure 1D).
10
11 Altogether, these data revealed that combination of MFC with BrdU incorporation or
12
13 KI67/DAPI staining is a practical method for MM cell proliferation assessment.
14
15 Furthermore, this methodology predicts EFS in newly diagnosed MM patients.
16
17

18 Discussion

19
20 During the last years, a significant expansion of MM therapeutic options occurred.
21
22 This underlines the need of useful biological and prognostic informations to develop
23
24 rational and targeted therapies. The plasma cell labeling index was shown to be one
25
26 of the most powerful risk factor in MM (10,11). However, this approach did not
27
28 separate malignant MM cells from normal PC. In the present study, we have
29
30 described that combination of multiparameter flow cytometry with BrdU incorporation
31
32 or ki67 staining allows specific MM cell proliferation assessment. We identified that
33
34 MM patients with $> 1.42\%$ proliferating MM cells identified with BrdU/DAPI or $> 1.1\%$
35
36 using ki67/DAPI correlated with significant quicker relapse, in the transplantation
37
38 settings. The methodology yields important information regarding myeloma prognosis
39
40 at the time of diagnosis. In this study, as expected, results obtained with BrdU
41
42 incorporation significantly correlated with ki67/DAPI analyses as described in other
43
44 cancers (27,28). Both approaches assess myeloma-specific proliferation on routinely
45
46 collected fresh bone marrow aspirates. Comparing the two methods, BrdU
47
48 incorporation requires ex vivo labeling of viable PC, in real time, with equipment that
49
50 could not be available in all clinical laboratories. The Ki67/DAPI methodology could
51
52 be used in any laboratories with equipment for clinical multiparametric flow cytometry
53
54
55
56
57
58
59
60

1
2
3 analyses. BrdU incorporation allows detection of S-phase cells only. Ki67 is found in
4 multiple phases of the cell cycle including S, G2, M (27,28) and could provide more
5 information relative to quiescence (29).
6
7

8
9 Several groups have shown the value of minimal residual disease (MRD) analysis
10 over previous response criteria in MM patients (3-5). Our methodology is also
11 suitable for MRD with low population of tumor cells. Combining the high clinical
12 sensitivity of MFC to the detection of a proliferative contingent of MM cells will enable
13 the identification of patients with active disease during the follow-up of MRD.
14
15

16 Furthermore, in the context of new effective myeloma treatment options, assessment
17 of MM cell proliferation may be valuable, in clinical trials, to identify novel agents that
18 could significantly affect the small proliferative compartment of MM cells.
19
20

21 In conclusion, combination of MFC with BrdU incorporation or ki67/DAPI staining
22 could become a standard for the determination of the fraction of proliferating
23 malignant plasma cells.
24
25
26
27
28
29
30
31
32
33
34
35

36 **Acknowledgements**

37 This work was supported by grants from French INCA (Institut National du Cancer)
38 Institute (2012-109/087437 and PLBIO15-256), Languedoc Roussillon CRLR
39 (R14026FF), LR-FEDER Hemodiag, ITMO Cancer (MM&TT) and AF3M.
40
41
42
43
44
45
46
47
48
49

50 **References**

- 51
52 1. Raab MS, Podar K, Breitkreutz I, Richardson PG, Anderson KC. Multiple myeloma.
53 Lancet 2009;374:324-39.
54 2. Rajkumar SV. Treatment of multiple myeloma. Nat Rev Clin Oncol 2011;8:479-91.
55
56
57
58
59
60

3. Landgren O, Owen RG. Better therapy requires better response evaluation: Paving the way for minimal residual disease testing for every myeloma patient. *Cytometry B Clin Cytom* 2016;90:14-20.
4. Paiva B, Almeida J, Perez-Andres M, Mateo G, Lopez A, Rasillo A, Vidriales MB, Lopez-Berges MC, Miguel JF, Orfao A. Utility of flow cytometry immunophenotyping in multiple myeloma and other clonal plasma cell-related disorders. *Cytometry B Clin Cytom* 2010;78:239-52.
5. Rawstron AC, Child JA, de Tute RM, Davies FE, Gregory WM, Bell SE, Szubert AJ, Navarro-Coy N, Drayson MT, Feyler S and others. Minimal residual disease assessed by multiparameter flow cytometry in multiple myeloma: impact on outcome in the Medical Research Council Myeloma IX Study. *J Clin Oncol* 2013;31:2540-7.
6. Maurer HR. Potential pitfalls of [3H]thymidine techniques to measure cell proliferation. *Cell Tissue Kinet* 1981;14:111-20.
7. Neckers LM, Funkhouser WK, Trepel JB, Cossman J, Gratzner HG. Significant non-S-phase DNA synthesis visualized by flow cytometry in activated and in malignant human lymphoid cells. *Exp Cell Res* 1985;156:429-38.
8. Falini B, Canino S, Sacchi S, Ciani C, Martelli MF, Gerdes J, Stein H, Pileri S, Gobbi M, Fagioli M and others. Immunocytochemical evaluation of the percentage of proliferating cells in pathological bone marrow and peripheral blood samples with the Ki-67 and anti-bromo-deoxyuridine monoclonal antibodies. *Br J Haematol* 1988;69:311-20.
9. San Miguel JF, Garcia-Sanz R, Gonzalez M, Moro MJ, Hernandez JM, Ortega F, Borrego D, Carnero M, Casanova F, Jimenez R and others. A new staging system for multiple myeloma based on the number of S-phase plasma cells. *Blood* 1995;85:448-55.
10. Greipp PR, Katzmann JA, O'Fallon WM, Kyle RA. Value of beta 2-microglobulin level and plasma cell labeling indices as prognostic factors in patients with newly diagnosed myeloma. *Blood* 1988;72:219-23.
11. Larsen JT, Chee CE, Lust JA, Greipp PR, Rajkumar SV. Reduction in plasma cell proliferation after initial therapy in newly diagnosed multiple myeloma measures treatment response and predicts improved survival. *Blood* 2011;118:2702-7.
12. Alexandrakis MG, Passam FH, Kyriakou DS, Dambaki K, Niniraki M, Stathopoulos E. Ki-67 proliferation index: correlation with prognostic parameters and outcome in multiple myeloma. *Am J Clin Oncol* 2004;27:8-13.
13. Alaterre E, Raimbault S, Garcia JM, Reme T, Requirand G, Klein B, Moreaux J. Automated and simplified identification of normal and abnormal plasma cells in Multiple Myeloma by flow cytometry. *Cytometry B Clin Cytom* 2017.
14. Hose D, Reme T, Meissner T, Moreaux J, Seckinger A, Lewis J, Benes V, Benner A, Hundemer M, Hielscher T and others. Inhibition of aurora kinases for tailored risk-adapted treatment of multiple myeloma. *Blood* 2009;113:4331-40.
15. Kassambara A, Reme T, Jourdan M, Fest T, Hose D, Tarte K, Klein B. GenomicScape: an easy-to-use web tool for gene expression data analysis. Application to investigate the molecular events in the differentiation of B cells into plasma cells. *PLoS Comput Biol* 2015;11:e1004077.
16. Gentleman RC, Carey VJ, Bates DM, Bolstad B, Dettling M, Dudoit S, Ellis B, Gautier L, Ge Y, Gentry J and others. Bioconductor: open software development for computational biology and bioinformatics. *Genome Biol* 2004;5:R80.

17. Subramanian A, Tamayo P, Mootha VK, Mukherjee S, Ebert BL, Gillette MA, Paulovich A, Pomeroy SL, Golub TR, Lander ES and others. Gene set enrichment analysis: a knowledge-based approach for interpreting genome-wide expression profiles. *Proc Natl Acad Sci U S A* 2005;102:15545-50.
18. Hothorn T, Lausen B. On the exact distribution of maximally selected rank statistics. *Computational Statistics & Data Analysis* 2003;43:121-137.
19. Bruyer A, Maes K, Herviou L, Kassambara A, Seckinger A, Cartron G, Reme T, Robert N, Requirand G, Boireau S and others. DNMTi/HDACi combined epigenetic targeted treatment induces reprogramming of myeloma cells in the direction of normal plasma cells. *Br J Cancer* 2018.
20. Kassambara A, Hose D, Moreaux J, Walker BA, Protopopov A, Reme T, Pellestor F, Pantesco V, Jauch A, Morgan G and others. Genes with a spike expression are clustered in chromosome (sub)bands and spike (sub)bands have a powerful prognostic value in patients with multiple myeloma. *Haematologica* 2012;97:622-30.
21. Viziteu E, Klein B, Basbous J, Lin YL, Hirtz C, Gourzones C, Tiers L, Bruyer A, Vincent L, Grandmougin C and others. RECQ1 helicase is involved in replication stress survival and drug resistance in multiple myeloma. *Leukemia* 2017;10:2104-2113.
22. Zhan F, Huang Y, Colla S, Stewart JP, Hanamura I, Gupta S, Epstein J, Yaccoby S, Sawyer J, Burington B and others. The molecular classification of multiple myeloma. *Blood* 2006;108:2020-8.
23. Hose D, Reme T, Hielscher T, Moreaux J, Messner T, Seckinger A, Benner A, Shaughnessy JD, Jr., Barlogie B, Zhou Y and others. Proliferation is a central independent prognostic factor and target for personalized and risk-adapted treatment in multiple myeloma. *Haematologica* 2011;96:87-95.
24. Shaughnessy JD, Jr., Zhan F, Burington BE, Huang Y, Colla S, Hanamura I, Stewart JP, Kordsmeier B, Randolph C, Williams DR and others. A validated gene expression model of high-risk multiple myeloma is defined by deregulated expression of genes mapping to chromosome 1. *Blood* 2007;109:2276-84.
25. Reme T, Hose D, Theillet C, Klein B. Modeling risk stratification in human cancer. *Bioinformatics* 2013;29:1149-57.
26. Decaux O, Lode L, Magrangeas F, Charbonnel C, Gouraud W, Jezequel P, Attal M, Harousseau JL, Moreau P, Bataille R and others. Prediction of survival in multiple myeloma based on gene expression profiles reveals cell cycle and chromosomal instability signatures in high-risk patients and hyperdiploid signatures in low-risk patients: a study of the Intergroupe Francophone du Myelome. *J Clin Oncol* 2008;26:4798-805.
27. Goodson WH, Moore DH, Waldman FM. Ki-67 correlates with in vivo bromodeoxyuridine labeling index in operable breast cancer. *J Clin Oncol* 2006;24:3809; author reply 3809.
28. Gasparri F, Wang N, Skog S, Galvani A, Eriksson S. Thymidine kinase 1 expression defines an activated G1 state of the cell cycle as revealed with site-specific antibodies and ArrayScan assays. *Eur J Cell Biol* 2009;88:779-85.
29. Yao G. Modelling mammalian cellular quiescence. *Interface Focus* 2014;4:20130074.

Figure legends:

1
2
3 **Figure 1. Assessment of malignant plasma cell proliferation using BrdU**
4 **incorporation**
5

6
7
8 (A) The percentage of malignant proliferative MM cells was analyzed using BrdU
9 incorporation and multiparameter flow cytometry in bone marrow samples from 44
10 patients with MGUS, 153 newly diagnosed patients and 69 patients at relapse.

11
12
13
14 (B) A percentage of proliferating MM cells > 1.42% could predict for shorter event
15 free survival. Patients of a cohort of 60 newly diagnosed patients were ranked
16 according to increasing proliferation defined by BrdU incorporation/DAPI and a
17 maximum difference in event free survival (EFS) was obtained using the Maxstat R
18 function.
19

20
21
22
23 (C) Top gene sets significantly associated with high risk proliferation assessment in
24 MM.

25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
GSEA enrichment plots with the absolute enrichment P value and the normalized
enrichment score of the gene set.

(D) Top gene sets significantly associated with low risk proliferation assessment in
MM.

GSEA enrichment plots with the absolute enrichment P value and the normalized
enrichment score of the gene set.

Figure 2. Assessment of malignant plasma cell proliferation using ki67/DAPI

(A) The percentage of proliferating malignant plasma cells was analyzed by flow
cytometry using ki67+/DAPI+ staining, in bone marrow samples from 11 patients with
MGUS, 49 newly diagnosed patients and 18 patients at relapse.

(B) A percentage of proliferating MM cells defined by ki67+/DAPI+ staining > 1.1%
could predict for shorter event free survival. Patients of a cohort of 31 newly
diagnosed patients were ranked according to increasing proliferation defined by

1
2
3 ki67+/DAPI+ staining and a maximum difference in event free survival (EFS) was
4
5 obtained using the Maxstat R algorithm.
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

For Peer Review

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60






Figure 1

Figure 1

150x299mm (300 x 300 DPI)



	Median survival (95% CI)
Ki67 ⁺ /DAPI ⁺ ≤ -1.1%	Not reached
Ki67 ⁺ /DAPI ⁺ > 1.1%	19 months (13.8; 24.1)

Figure 2

Figure 2

194x274mm (300 x 300 DPI)



1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47

Tubes	DAPI	FITC	PE	ECD	PeCy5.5	PeCy7	APC	APC-AF750
1 – Control	DAPI	Lambda	Kappa	CD45	CD38	CD19	Anti BrdU	CD20
2 - Negative Pool	DAPI	Lambda	Kappa	CD45	CD38	Negative Pool CD19/CD27	Anti BrdU	CD20
2 - Positive Pool	DAPI	Lambda	Kappa	CD45	CD38	Positive Pool CD56/CD117/CD200	Anti BrdU	CD20


Table 1A: BrdU panel composition. The panel comprises 3 tubes that contain DAPI, anti-CD45, anti-CD20, anti-CD38, anti-BrdU, anti-kappa, anti-lambda antibodies, a positive antibody pool against markers that are overexpressed in malignant PC (CD56, CD117 and CD200) and a negative antibody pool against markers that are downregulated in malignant PC (CD19 and CD27) compared to normal PC.

Tubes	DAPI	FITC	PE	ECD	PeCy5.5	PeCy7	APC	APC-AF750
1 – Control	DAPI	Lambda	Kappa	CD45	CD38	CD19	Ki67	CD20
2 - Negative Pool	DAPI	Lambda	Kappa	CD45	CD38	Negative Pool CD19/CD27	Ki67	CD20
2 - Positive Pool	DAPI	Lambda	Kappa	CD45	CD38	Positive Pool CD56/CD117/CD200	Ki67	CD20

Table 1B: Ki67 panel composition. The panel comprises 3 tubes that contain DAPI, anti-CD45, anti-CD20, anti-CD38, anti-Ki67, anti-kappa, anti-lambda antibodies, a positive antibody pool against markers that are overexpressed in malignant PC (CD56, CD117 and CD200) and a negative antibody pool against markers that are downregulated in malignant PC (CD19 and CD27) compared to normal PC.



Supplementary Figure S1A



Supplementary Figure S1B


1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47

Supplementary Figure S1. Immunophenotyping of normal and malignant plasma cells.

(A) Singlets are selected on FSC and SSC plots and total leukocytes on FSC/SSC plots. PC and B cells are selected on CD45/CD38 and CD20/CD19 plots. Events on the kappa/lambda diagonal are removed. Abnormal PC are selected on the CD27, CD56, CD117, and CD200 dimensions (Alaterre, et al. 2017).



(B) Assessment of proliferating MM cells using BrdU incorporation and ki67 immunostaining.

Alaterre E, Raimbault S, Garcia JM, Reme T, Requirand G, Klein B, Moreaux J. Automated and simplified identification of normal and abnormal plasma cells in Multiple Myeloma by flow cytometry. *Cytometry B Clin Cytom* 2017.




Supplementary Figure S2: Correlation between combined MFC/BrdU incorporation and MFC/Ki67+/DAPI+ staining for malignant plasma cell proliferation assessment. Proliferation has been investigated by flow cytometry in 11 MGUS and 49 MM patients at diagnosis.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47




Supplementary Figure S3: Proliferation staining combining (A) MFC and BrdU incorporation (n=81) or (B) MFC and Ki67 staining (n=30) in the 7 subgroups of the UAMS molecular classification of multiple myeloma. Proliferation staining was investigated in newly-diagnosed MM patients. PR: proliferation, LB: low bone disease, MS: MMSET, HY: hyperdiploid, CD1: Cyclin D1-Cyclin D3, CD2: Cyclin D1-Cyclin D3, MF: MAF, MY: myeloid. * Indicate that proliferation staining is significantly higher in the group compared to all the patients of the cohort ($P < 0.05$).


A




B



C



D



Supplementary Figure S4.

The percentage of proliferating MM cells using BrdU was assessed for patients classified according to GEP-based GPI index (Hose, *et al* 2011) (A), HRS (Shaughnessy, *et al* 2007) (B), RS (Reme, *et al* 2013) (C) and IFM (Decaux, *et al* 2008) (D) risk scores (n = 83). * indicates a P-value <0.05 using a Student t-test.


Hose, D., Reme, T., Hielscher, T., Moreaux, J., Messner, T., Seckinger, A., Benner, A., Shaughnessy, J.D., Jr., Barlogie, B., Zhou, Y., Hillengass, J., Bertsch, U., Neben, K., Mohler, T., Rossi, J.F., Jauch, A., Klein, B. & Goldschmidt, H. (2011) Proliferation is a central independent prognostic factor and target for personalized and risk-adapted treatment in multiple myeloma. *Haematologica*, **96**, 87-95.

Shaughnessy, J.D., Jr., Zhan, F., Burington, B.E., Huang, Y., Colla, S., Hanamura, I., Stewart, J.P., Kordsmeier, B., Randolph, C., Williams, D.R., Xiao, Y., Xu, H., Epstein, J., Anaissie, E., Krishna, S.G., Cottler-Fox, M., Hollmig, K., Mohiuddin, A., Pineda-Roman, M., Tricot, G., van Rhee, F., Sawyer, J., Alsayed, Y., Walker, R., Zangari, M., Crowley, J. & Barlogie, B. (2007) A validated gene expression model of high-risk multiple myeloma is defined by deregulated expression of genes mapping to chromosome 1. *Blood*, **109**, 2276-2284.

Reme, T., Hose, D., Theillet, C. & Klein, B. (2013) Modeling risk stratification in human cancer. *Bioinformatics*, **29**, 1149-1157.


Decaux, O., Lode, L., Magrangeas, F., Charbonnel, C., Gouraud, W., Jezequel, P., Attal, M., Harousseau, J.L., Moreau, P., Bataille, R., Campion, L., Avet-Loiseau, H. & Minvielle, S. (2008) Prediction of survival in multiple myeloma based on gene expression profiles reveals cell cycle and chromosomal instability signatures in high-risk patients and hyperdiploid signatures in low-risk patients: a study of the Intergroupe Francophone du Myelome. *J Clin Oncol*, **26**, 4798-4805

A




1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47

B




C



1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47

D



Supplementary Figure S5.

The percentage of proliferating MM cells was assessed for patients classified according to GEP-based GPI index (A), HRS (B), RS (C) and IFM (D) risk scores (n = 31). * indicates a P-value <0.05 using a Student t-test.